Feline Pancreatic Lipase Test Kit  
For veterinary use only.

English version

Feline SNAP® fPL\*  
In vitro test for the determination of pancreas-specific lipase levels in feline serum.

**Precautions and warnings**

- All wastes should be properly decontaminated prior to disposal.
- Do not mix components from kits with different lot numbers.
- Do not use a SNAP device that has been activated prior to the addition of sample.
- Refer to country specific Material Safety Data Sheet for regional hazard identification.

**Storage**

- Store at 2-8°C until expiration date.

**Kit components**

Item	Reagents	Quantity
1	1 bottle of anti-chicken:HRPO/anti-fPL-HRPO conjugate	4,0 mL
2	SNAP device	5 or 10

Each SNAP device contains 0.4 mL of wash solution and 0.6 mL of substrate solution.

**Other components:** transfer pipettes, sample tubes, and reagent rack

**Conjugate** - May cause an allergic skin reaction. Harmful to aquatic life with long lasting effects. Avoid breathing mist/vapors. Wear protective gloves. If skin irritation or rash occurs: Get medical advice/attention.

**Wash Solution** - May cause an allergic skin reaction. Avoid breathing mist/vapors. Wear protective gloves. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention.

**Sample information**

- Samples must be at room temperature (15-25°C) before beginning the test procedure.
- Serum, either fresh or stored at 2-8°C for up to one week, can be used.
- For longer storage, serum can be frozen (-20°C or colder) and then re centrifuged before use.
- Hemolyzed, icteric or lipemic samples will not affect test results.

**Test procedure**

- Allow all components to equilibrate at room temperature (15-25°C) for at least 30 minutes. **Do not heat.**

- Hold the required transfer pipette vertically and dispense 3 drops of sample into a new sample tube.

- Holding the bottle vertically, add 4 drops of conjugate to the sample tube.

- Cap the sample tube and mix thoroughly by inverting it 3-5 times.

- Place the device on a horizontal surface. Add the entire contents of the sample tube to the sample well, being careful not to splash the contents outside of the sample well.

The sample will flow across the result window, reaching the activation circle in 30-60 seconds. Some sample may remain in the sample well.


- When color **FIRST** appears in the activation circle, push the activator firmly until it is flush with the device body.


**NOTE:** Some samples may not flow to the activation circle within 60 seconds. In this case, press the activator after the sample has flowed across the result window.

- Read test result **10 minutes** from the time of activation.

**Interpreting the test results**

**Abnormal level**  
If the color intensity of the sample spot is equal to or darker than the color intensity of the reference spot, fPL levels are abnormal.

**Normal level**  
If the color intensity of the sample spot is lighter than the color intensity of the reference spot, fPL levels are normal.  
**NOTE:** It is possible for no color to develop on the sample spot.

**Invalid results**


- Background**—If the sample is allowed to flow past the activation circle, background color may result. Some background color is normal. However, if colored background obscures test result, repeat the test.
- No color development**—If reference spot does not develop color, repeat the test.

**IDEXX SNAPshot Dx\* Analyzer**


Test results can also be read using the SNAPshot Dx\* Analyzer. A complete description of how to enter patient data and read test results using the SNAPshot Dx Analyzer can be found in the SNAPshot Dx Analyzer user guide.

Symbol descriptions	Use by date	EC/REP	Authorized representative in the European Community
	Batch code (lot)		
	Serial number		Consult instructions for use
	Temperature limitation		In vitro diagnostic
	Manufacturer		Date of manufacture
	Catalog number		

IDEXX Technical Support  
USA/Canada: 1 800 248 2483 • idexx.com.au  
Australia: 1300 44 33 99 • idexx.com.au  
Europe: idexx.eu

\* SNAP and SNAP fPL are trademarks or registered trademarks of IDEXX Laboratories, Inc. or its affiliates in the United States and/or other countries. One IDEXX Drive Westbrook, Maine 04092 USA idexx.com © 2022 IDEXX Laboratories, Inc. All rights reserved. • 06-2852-06

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