




ImageVue Peripherals and Accessories

Features	Benefits	Dimensions	What's in the box?	
System computer—all ImageVue* systems				
 <p data-bbox="159 618 417 643">Dell* OptiPlex* XE3 SFF PC</p>	<ul style="list-style-type: none"> • Processor: 8th Gen Intel* Core* i5-8500 • RAM: 16 GB • Hard drive: 1 TB (7,200 rpm) • Windows* 10 Pro 64-bit operating system 	<ul style="list-style-type: none"> • Ready to connect to in-practice networks or the web 	<ul style="list-style-type: none"> • Height: 11.42" • Width: 3.7" • Depth: 11.5" • Weight: 11.57 lb 	<ul style="list-style-type: none"> • Computer • Standard wired mouse • Standard wired keyboard • Power cable
Monitor—ImageVue digital radiography (DR) systems				
 <p data-bbox="121 915 459 938">Elo* 2201L 22" Touchscreen Monitor</p>	<ul style="list-style-type: none"> • Aspect ratio: 16:9 • Display area: 18.5" x 10.63" • Resolution: 1920 x 1080 @ 60 Hz 	<ul style="list-style-type: none"> • Touchscreen capability • Wide viewing angle 	<ul style="list-style-type: none"> • Height: 21.7" • Width: 15.4" • Depth: 7.6" • Weight: 15.0 lb 	<ul style="list-style-type: none"> • Touchscreen monitor • Stand (riser and base) • Power cable • Cables • Quick install guide
Monitor—ImageVue computed radiography (CR) system				
 <p data-bbox="100 1240 478 1284">Barco* MDRC-1219 TS 19" Touchscreen Monitor</p>	<ul style="list-style-type: none"> • Aspect ratio: 5:4 • Display area: 14.8" x 11.8" • Resolution: 1280 x 1024 @ 60 Hz 	<ul style="list-style-type: none"> • Touchscreen capability • Compact monitor footprint 	<ul style="list-style-type: none"> • Height: 13.7" • Width: 17.3" • Depth: 7.9" (with stand) • Weight: 13.4 lb 	<ul style="list-style-type: none"> • Touchscreen monitor • Stand (riser and base) • Power cable • Cables • Quick setup guide